


Registered

KF94  at KFDA

CARE ALL

Yellow Dust · Disinfection Mask

High Protection
rate against dust
particles, Hazardous
Material and sources
of the respiratory
infection around
over 94%



Outside and inside filters
exclusively provided
Yu-Han Kimberly Co., Ltd
who is the best manufacturer in Korea.
All of other material is made in Korea,
like nose fixed frame,



Dust particles, Hazardous Material and sources of the respiratory infection must be protected by **outside & inside Filter** which are made with high quality and effectiveness by

Yu-Han Kimberly Co., Ltd. in Mark

Yu-Han Kimberly had a named Company who could manufacture baby diapers, it perfectly remove the skin allergy and troubles. Apple health Care Co., Ltd (AHC) can serve mask with ergonomic technics and high-performed comfortable wearing.

AHC use the middle of Mask which is the core of mask made high-density & **ultra-static MB filters**. The one of AHC's MB filter passed Dust Collection Efficiency by **paraffin oil test** (highest level-covered respiratory sources) in KFDA (over **94%** rates). It proves AHC's Mask performance against dust particles, Hazardous Material and sources of the respiratory infection.

This nose wire were made by non-hazardous material which maximizes the **cohesive** forces to serve excellent **wearables**.

It causes to protect the zero-tolerance of respiratory leakages and the effective **antifogging** for whom somebody had a glasses.

For reason that ear string is the key product point because ear string covers the resilience force when mask's wearing period, AHC use only Korea Origin and the material for baby which had been hardly passed for by Korean Safety Regulation. **Specially ear string, which is two items(2) of pH, pormaldehyde, Twenty-four(24) items of Nonylphenol, Twenty-three(23) items of allergy caused dyes, Total Fifty-one(51) items of hazardous material** had been tested and certified by Korea intertek who is international certifier for all chemical material. It definitely serves non-allergy pure and comfortable wearing to customers.

“Made in Korea - 100%”

This is AHC's Sprit for the Product.
Did you check ALL of material which is

“Made in Korea - 100%”?

Did think that you will buy masks which made by low-qualified material and no efficient rate one?



Why some masks has 4-layers or 5-layers?
AHC ensures that care-All is the
Mask by KFDA's certified KF94,
3-layers, Made-in-Korea one.

AHC cannot does not want to compare to the products which made heavy, low-material, 4 or 5 layers masks. AHC NERVE accepts it. Due to too much demands for the public health purposes, there are a lots of products sold out with low prices in public, which made by chinese's low quality MB filter mask. Every public did not acknowledged well, this MB filter only can protect respiratory virus.



care-all is 100% three layers made in Korea Mask.
AHC is very pride of manufacturing the art of technology, ergonomic design, and wearing easy product.

“Please wear AHC's Mask”

Experience the difference first handed.

Now, you are wearing the mask which made that is china origin material, made in Korea?



Certification by KFDA

As nonmedical product which is licensed
product passed by complicated & strict test

✓ Manufacturing Report Card



✓ Item License



✓ Test Report

Test Result for Dust Collection Efficiency Test by Sodium Chloride (98%/99%/99%)
Test Result for Dust Collection Efficiency Test by Paraffin Oil (98%/99%/99%)



✓ Test Report for ear string



✓ Certification by FDA



How to use the product



- 1 Check that the inside of the mask is not contaminated.
- 2 Wear the mask to cover your nose and mouth.
- 3 Hook each strap over your ear to lock them in place.
- 4 Cover the entire mask with both hands and adjust it to fit the face to make sure there is no air leakage.

Precautions for use

- Make sure not to place towel, tissue, etc. inside of the mask.
- Do not use if the inside of the mask is contaminated.
- Do not wash for use.
- Do not crush or deform the mask.
- Do not touch the surface of the mask after wearing.
- If you are a pregnant women, a child, a senior citizen, a patient suffering from respiratory tract or cardiovascular disease and feel uncomfortable in breathing when you wear the mask please discontinue using immediately and consult your doctor or other specialist if necessary.

Quality Indicators

Product Name	CARE ALL KF94		
Ingredients	nonwoven(outer, filter, lining), plastic coated wire, nylon string		
Manufacturer	APPLE HEALTH CARE Co., Ltd. / 53, Myeongbonggeonam-ro, Onsan-eup, Ulju-gun, Ulsan, Republic of Korea		
Distributor	APPLE HEALTH CARE Co., Ltd. / 53, Myeongbonggeonam-ro, Onsan-eup, Ulju-gun, Ulsan, Republic of Korea		
	DB HOLDINGS Co., Ltd. 18-1, Samsan-ro 392beon-gil, Nam-gu, Ulsan, Republic of Korea		
Country of manufacture	KOREA	Type of Mask	Disinfection Mask
Contact Us	Tel. +82.52.239.2016 / Fax. +82.52.239.2017		

If there is something wrong with this product, it can be exchanged and refunded according to the consumer dispute resolution standard announced by the Fair Trade Commission.



APPLE HEALTH CARE